

AUTO CHART PROJECTOR (LED LAMP)

Features:

- 1.By the way of double locating, it would correct errors automatically.
- 2.Offering more clearer and sharper pictures.
- 3.Increase the picture density; there are 5 icons in one line.
- 4.The charts are bigger,size:33*22.5(CM)(in 5M refraction).
- 5.Projected by wireless remote controller.
- 6.Able to project the charts one by one.
- 7.Standard accessories, the metal board is covered by bias material.



Specifications:

- 1.Projection distance:2-6M
- 2.Projection magnification:30×(in 5M refraction)
- 3.Projection size :33*22.5(CM),Φ275MM(in 5M refraction)
- 4.Numbers of charts:33
- 5.Chart change-over:3 frame/sec
- 6.Numbers of masks:1 open,3 horizontal line,8 vertical line, 16 single isolation,1 red and green
- 7.Tilt range:10°horizontally as well as vertically
- 8.Bulb:high-luminance and long-life LED light source
- 9.Image shift interval:0.5 second
- 10.Automatic shut-off mechanism:5 minutes
- 11.Input voltage: AC 110V/60HZ or 220V/50HZ
- 12.Power consumption:65W
- 13.Dimensions :19.6(W)*29.2(L)*23.4(H)(CM)
- 14.Weight:6.5KGS

Charts and alternative accessories:

E 4.18	MEW 4.18 3ME 4.18	MEWME 4.18 3EME 4.18 WME3 4.18	EWWE 4.18 W3E 4.18 W3EM 4.18	W3WME 4.18 EM3W 4.18 M3E3M 4.18	W3WME 4.18 EM3W 4.18 M3E3M 4.18	FVN 4.18 ZPD 4.18
TDEPA 4.18 VLNEK 4.18 DAOFR 4.18	ECKHD 4.18 FZDKV 4.18 OFLCN 4.18	APEDT 4.18 TZVEA 4.18 OHPNF 4.18	KENLV 4.18 EVOTD 4.18 VKDZF 4.18	8 4.18	592 4.18 638 4.18	35698 4.18 68239 4.18 86953 4.18
23586 4.18 69328 4.18 92563 4.18	38965 4.18 59682 4.18 95836 4.18	63589 4.18 86395 4.18 39826 4.18	4.18	4.18	4.18	4.18
4.18	4.18	53986 35869	4.18	4.18	4.18	4.18
4.18	4.18	4.18	4.18	4.18	4.18	CHARTS



CT1502
Wall mount bracket



CT1503



CT1503A

Stand for auto chart projector